



No. : WZH054701005K-2

TEST REPORT

Sample Name: TILT DISPCAY LABEL

Applicant Name: Wuxi Prmierse Technology Co,. Ltd

Manufacturer: Wuxi Prmierse Technology Co,. Ltd

Supplier: /

Buyer: /

Shenzhen Yitong Testing Technology Co., Ltd





TEST REPORT

No.: WZH054701005K-2

Date: Nov 21, 2022

Page 1 Of 4

Applicant Name: Wuxi Prmierse Technology Co,. Ltd

Address: No.210 Xinzhou Road Meicun town Xinwu District Wuxi City

Sample Information

Sample Name: TILT DISPCAY LABEL

Sample No.: WZH054701005-1 to 2

Sample Model: TILTSEE

Sample Size: 74 mm×59 mm×5 mm

Sample Weight: 0.0122 kg

Sample Quantity: 2 PCS

Sample Description: Normal

Sample Supplied By: Client

Manufacturer: Wuxi Prmierse Technology Co,. Ltd

Test Information

Test Items: Low Temperature Test, High Temperature Test

Receive Date: Nov 14, 2022

Testing Date: Nov 14, 2022 to Nov 15, 2022

Preconditioning: /

Test Result

Sample No.	Test Items	Test Result
WZH054701005-1	Low Temperature Test	After the test, there is no appearance abnormality on the sample by visual inspection.
WZH054701005-2	High Temperature Test	After the test, there is no appearance abnormality on the sample by visual inspection.

Authorized: Dicky



No.: WZH054701005K-2

Date: Nov 21, 2022

Page 2 Of 4

Test Item 1: Low Temperature Test

1. Test Standard: Supplied by client

2. Test Equipment:

Equipment Name	Equipment Model	Equipment No.	Calibration Useful-life
Temperature & Humidity Chamber	GDJS-500D	TTS-YQ-197	Apr 10, 2023

3. Test Remark & Notes:

- (1) Place the sample into the Temperature & Humidity Chamber;
- (2) Adjust the inner temperature to -10 °C, keep the sample in that environment for 24 h;
- (3) Check the sample after the test.

Test Item 2: High Temperature Test

1. Test Standard: Supplied by client

2. Test Equipment:

Equipment Name	Equipment Model	Equipment No.	Calibration Useful-life
Temperature & Humidity Chamber	GDJS-408C	TTS-YQ-196	Apr 10, 2023

3. Test Remark & Notes:

- (1) Place the sample into the Temperature & Humidity Chamber;
- (2) Adjust the inner temperature to 60 °C, keep the sample in that environment for 24 h;
- (3) Check the sample after the test.

Appendix: Test Pictures

Pre-test:



sample WZH054701005-1




sample WZH054701005-2

Testing:

Low Temperature Test:



High Temperature Test:



Post-test:

sample WZH054701005-1:



No appearance abnormality on the sample by visual inspection

sample WZH054701005-2:



No appearance abnormality on the sample by visual inspection

- - - END OF THE REPORT - - -

The test report with no signature or specialized stamp is invalid. The result(s) shown in this report refer only to the sample(s) tested. Without the approval of TTS, this report can't be reproduced except in full. The testing date and result(s) in this report is(are) just for scientific research, education, enterprise product development and internal quality control etc.